


Product Profile

Part number	1808-1022
Article name	Adapter for microplate . Flat bottom High Base
Category	Standard
CO Customs tariff code	DE 8479 90 70

Specifications

Ideal fit for*	1536 SW PS F-bottom, HiBase, solid/ μ Clear/glass bo Greiner Bio One 782061, 782073, 782074, 782075, 782076, 782077, 782078, 782080, 782086, 782092, 782093, 782094 384 SW PP Small Volume TM HiBase, solid bottom Greiner Bio One 784075, 784076, 784101, 784075-25
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 1808-xxx (All instruments with three point connection)
Mounting material	3x Cylinder screw (BN20146 M3x6-A2, lower head) 3x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	165 g 5.82 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!