


**Product Profile**

Part number	2016-1022
Article name	Microplate adapter - Flat bottom, high base
Category	Standard
CO   Customs tariff code	DE   8479 90 70

**Specifications**

Ideal fit for*	1536-well Black Polystyrene NBS Microplate PS Corning   3724, 3725, 3726, 3727, 3728, 3729  1536 SW PS F-bottom, HiBase, solid/ $\mu$ Clear/glass bo Greiner Bio One   782061, 782073, 782074, 782075, 782076, 782077, 782078, 782080, 782086, 782092, 782093, 782094  384 SW PP Small Volume <sup>TM</sup> HiBase, solid bottom Greiner Bio One   784075, 784076, 784101, 784075-25
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	122 g   4.3 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!