



Product Profile

Part number	2016-1041
Article name	Adapter for PCR Plate . 96 well
Category	PCR microplates
CO Customs tariff code	DE 8479 90 70

Specifications

Ideal fit for*	<p>PCR plate, 96-well, thin-wall, skirted, low profil Ritter 43001-0302</p> <p>Axygen® 96-well PP PCR Microplate, no skirt, clear Corning </p> <p>twin.tec® PCR Plate 96 semi-skirted 250µl Eppendorf 0030128575, 0030128583, 0030128591, 0030128605, 0030128613, 0030129326, 0030129334, 0030129610</p>
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	206 g 7.27 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!