



**Product Profile**

Part number	2016-1181
Article name	Adapter - Sarstedt® MegaBlock 96 Well 2.2 ml
Category	Deep well & storage plates
CO   Customs tariff code	DE   8479 90 70

**Specifications**

Ideal fit for*	Sarstedt Deep Well MegaBlock®, 0,5 ml, PP Sarstedt   Deep Well MegaBlock®, 96 Well, 1,2 ml, PS, Rundbod Sarstedt   82.1971.002 in PP Deep Well MegaBlock®, 96 Well, 2,2 ml, PP, Rundbod Sarstedt
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	Qinstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	114 g   4.02 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!