


Product Profile

Part number	2016-1217
Article name	Adapter for 96/450µl microplate
Category	Standard
CO Customs tariff code	DE 8479 90 70

Specifications

Ideal fit for*	Riplate® 96 SRW magnetic, 2 ml Ritter 96 DW PP MASTERBLOCK®, 1 ml Greiner Bio One 780201, 780203, 780204, 780205, 780206, 780215, 780261, 780263, 780264, 780265, 780266 96 Square Deep Well Microplate, KingFisher™ Style Azenta
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	Qinstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	124 g 4.37 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!

