

Slide Rack

Product Information

The **Slide Rack** is a special rack for slides with the standard size of 1" x 3". It is a ideal tool where to place your up to 8 slides!

Advantages

- Convenient **parallel handling** for up to **8 slides**
- Ideal storage for assembling and sealing of slides
- **Prevents scratching** of the bottom
- Ideal storage between and after processing
- Allows **easy gas exchange** through the bottom of the slides
- **Stackable**
- **Made of aluminium anodised**

Technical specifications

Dimensions over all	210 x 85 mm
Number of slides	1 - 8
Slide dimensions	Supports multiple platforms (1" x 3") Length: 75.0 - 77.5 mm Width: 25.0 - 26.0 mm
Material	All-aluminium anodised

Order Information

Name	Article number
Slide Rack	1301-0011

Your contact

QUANTIFOIL Instruments GmbH
Loebstedter Str. 101
07749 Jena
Germany

Phone: +49 (0)3641-87612-0
Fax: +49 (0)3641-87612-99
www.QInstruments.com
info@QInstruments.com



Slide Rack



**Slide Rack
with different
types of slides**

