


Product Profile

Part number	1808-1032
Article name	Adapter for microplate . 96 well round bottom
Category	Standard
CO Customs tariff code	DE 8479 90 70

Specifications

Ideal fit for*	StorPlate-96 U, DW, (2.0mL), DW PP, 96 (U), 2.0ml Perkin Elmer Zellkulturplatte, 96 Well, Rundboden, PS Sarstedt 83.3925, 83.3925.400, 83.3925.500, Riplate® medio, PP, 2ml Ritter
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	Qinstruments Devices with Art No 1808-xxx (All instruments with three point connection)
Mounting material	3x Cylinder screw (BN20146 M3x6-A2, lower head) 3x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	127 g 4.48 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!

