



**Product Profile**

Part number	1808-1041
Article name	Adapter for PCR Plate . 96 well
Category	PCR microplates
CO   Customs tariff code	DE   8479 90 70

**Specifications**

Ideal fit for*	PCR plate, 96-well, thin-wall, skirted, low profil Ritter   43001-0302  Axygen® 96-well PP PCR Microplate, no skirt, clear Corning
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 1808-xxx (All instruments with three point connection)
Mounting material	3x Cylinder screw (BN20146 M3x6-A2, lower head)   3x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	228 g   8.04 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!