

Product Profile


Part number	2016-1021
Article name	Microplate adapter - Flat bottom standard
Category	Standard
CO Customs tariff code	DE 8419 89 98

Specifications

Ideal fit for*	96 SW PS F-bottom / standard, solid bottom Greiner Bio One 655001, 655061, 655101, 655161, 651199 MICROPLATTE, 384 WELL, PS, F-BODEN Greiner Bio One 781073, 781074, 781075, 781076, 781077, 781079, 781080, 781086 sciPLEXPLATE 96 BICO CPH-5551-10 black Greiner 96
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	63 g 2.24 oz
Recommended max rpm**	5000

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!