


**Product Profile**

Part number	2016-1021
Article name	Microplate adapter - Flat bottom standard
Category	Standard
CO   Customs tariff code	DE   8479 90 70

**Specifications**

<p><b>Ideal fit for*</b></p> <p>96 SW PS F-bottom / standard, solid bottom Greiner Bio One   655001, 655061, 655101, 655161, 651199</p> <p>MICROPLATTE, 384 WELL, PS, F-BODEN Greiner Bio One   781073, 781074, 781075, 781076, 781077, 781079, 781080, 781086</p> <p>sciPLEXPLATE 96 BICO   CPH-5551-10 black Greiner 96</p>	<p><b>Material / Color</b></p> <p>Aluminum anodized</p>
<p><b>Temperature range</b></p> <p>-20 to 100 °C</p>	
<p><b>Compatible with</b></p> <p>QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)</p>	
<p><b>Mounting material</b></p> <p>1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8</p>	
<p><b>Adapter dimensions</b></p> <p>see drawing</p>	
<p><b>Weight</b></p> <p>64 g   2.26 oz</p>	
<p><b>Recommended max rpm</b></p> <p>The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.</p>	

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!