

Product Profile


Part number	2016-1022
Article name	Microplate adapter - Flat bottom, high base
Category	Standard
CO Customs tariff code	DE 8419 89 98

Specifications

<p>Ideal fit for*</p> <p>1536-well Black Polystyrene NBS Microplate PS Corning 3724, 3725, 3726, 3727, 3728, 3729</p> <p>1536 SW PS F-bottom, HiBase, solid/μClear/glass bo Greiner Bio One 782061, 782073, 782074, 782075, 782076, 782077, 782078, 782080, 782086, 782092, 782093, 782094</p> <p>384 SW PP Small VolumeTM HiBase, solid bottom Greiner Bio One 784075, 784076, 784101, 784075-25</p>	
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	122 g 4.3 oz
Recommended max rpm**	5000

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!