

## Adapter specification

---

**Product no.:** 2016-1022

**Product name:** Microplate adapter – Flat bottom hibase

---

**Description:** Adapter for cooling and heating of microplates  
Temperature range -10 to 100 °C  
Compatible with all BioShake/ColdPlate/HeatPlate automation modules  
with single point connection

**Ideal fit for:** Greiner 65590 and similar microplates

**Material / Color:** Aluminum anodized


**Only using for:** QInstruments Modules with Art. No. 2016-xxx

---

**Manufacturer name:** Quantifoil Instruments GmbH

**Address:** Loebstedter Str. 101 . 07749 Jena . Germany

---



**Drawing:**

