

Product Profile


Part number	2016-1032
Article name	Adapter for microplate . 96 well round bottom
Category	Standard
CO Customs tariff code	DE 8419 89 98

Specifications

<p>Ideal fit for*</p> <p>Material / Color</p> <p>Temperature range</p> <p>Compatible with</p> <p>Mounting material</p> <p>Adapter dimensions</p> <p>Weight</p> <p>Recommended max rpm**</p>	<p>Deep Well MegaBlock®, 96 Well, 1,2 ml, PS, Rundbod Sarstedt 82.1971.002 in PP</p> <p>Mikrotestplatte, 96 Well, Rundboden, PS Sarstedt 82.1582, 82.1582.100, 82.1582.200, 82.1582.001</p> <p>Deep Well MegaBlock®, 96 Well, 2,2 ml, PP, Rundbod Sarstedt </p> <p>Aluminum anodized</p> <p>-20 to 100 °C</p> <p>Qinstruments Devices with Art No 2016-xxx (All instruments with single point connection)</p> <p>1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8</p> <p>see drawing</p> <p>95 g 3.35 oz</p> <p>3000</p>
---	--

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!