

Product Profile



Part number	2016-1041
Article name	Adapter for PCR Plate . 96 well
Category	PCR microplates
CO Customs tariff code	DE 8419 89 98

Specifications

<p>Ideal fit for*</p> <p>Material / Color</p> <p>Temperature range</p> <p>Compatible with</p> <p>Mounting material</p> <p>Adapter dimensions</p> <p>Weight</p> <p>Recommended max rpm**</p>	<p>PCR plate, 96-well, thin-wall, skirted, low profil Ritter 43001-0302</p> <p>Axygen® 96-well PP PCR Microplate, no skirt, clear Corning </p> <p>twin.tec® PCR Plate 96 semi-skirted 250µl Eppendorf 0030128575, 0030128583, 0030128591, 0030128605, 0030128613, 0030129326, 0030129334, 0030129610</p> <p>Aluminum anodized</p> <p>-20 to 100 °C</p> <p>QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)</p> <p>1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8</p> <p>see drawing</p> <p>206 g 7.27 oz</p> <p>2500</p>
---	--

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!