

## Adapter specification

---

**Product no.:** 2016-1041

**Product name:** Adapter | 96PCR

---

**Description:** Adapter for cooling and heating of microplates  
Temperature range -10 to 100 °C  
Compatible with all BioShake/ColdPlate/HeatPlate automation modules  
with single point connection

**Ideal fit for:** Eppendorf® twin.tec® PCR Plate 96 (#0030 128.672) and similar plates

**Material / Color:** Aluminum anodized



**Only using for:** QInstruments Modules with Art. No. 2016-xxx

---

**Manufacturer name:** Quantifoil Instruments GmbH

**Address:** Loebstedter Str. 101 . 07749 Jena . Germany

---



**Drawing:**

