

**Product Profile**



Part number	2016-1042
Article name	Microplate adapter - 96 well Bio-Rad PCR plate
Category	PCR microplates
CO   Customs tariff code	DE   8419 89 98

**Specifications**

Ideal fit for*	PCR Platte, Vollrand, 96 Well, PP, 100µl Sarstedt   72.1980, 72.1980.010, 72.1980.201, 72.1980.600, 72.1980.700
	PCR-Platte, Halbrand, 96 Well, 200µl, PP, High Pro Sarstedt   72.1979, 72.1979.201, 72.1979.010, 72.1979.003, 72.1979.700, 72.1979.102, 72.1979, 132
	PCR-Platte Halbrand, 96 Well, transparent, 100µl Sarstedt   72.1981, 72.1981.010, 72.1981.202, 72.1981.232
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	223 g   7.87 oz
Recommended max rpm**	2500

\* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

\*\* The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!