



Product Profile

Part number	2016-1042
Article name	Microplate adapter - 96 well Bio-Rad PCR plate
Category	PCR microplates
CO Customs tariff code	DE 8479 90 70

Specifications

Ideal fit for*	PCR Platte, Vollrand, 96 Well, PP, 100µl Sarstedt 72.1980, 72.1980.010, 72.1980.201, 72.1980.600, 72.1980.700
	PCR-Platte, Halbrand, 96 Well, 200µl, PP, High Pro Sarstedt 72.1979, 72.1979.201, 72.1979.010, 72.1979.003, 72.1979.700, 72.1979.102, 72.1979, 132
	PCR-Platte Halbrand, 96 Well, transparent, 100µl Sarstedt 72.1981, 72.1981.010, 72.1981.202, 72.1981.232
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	223 g 7.87 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!