

## Adapter specification

---

**Product no.:** 2016-1042

**Product name:** Microplate adapter - 96 well Bio-Rad PCR plate

---

**Description:** Adapter for cooling and heating of microplates  
Temperature range -10 to 100 °C  
Compatible with all BioShake/ColdPlate/HeatPlate automation modules  
with single point connection

**Ideal fit for:** Bio-RAD (#MSP-9601, #MSP-9901) and similar disposables  
Fortitude #4TI-0960 (FrameStar® 96 Well Skirted PCR plate 200 µl)

**Material / Color:** Aluminum anodized




**Only using for:** QInstruments Modules with Art. No. 2016-xxx

---

**Manufacturer name:** Quantifoil Instruments GmbH

**Address:** Loebstedter Str. 101 . 07749 Jena . Germany

---



**Drawing:**

