


Product Profile

Part number	2016-1051
Article name	Adapter for PCR Plate . 384 well
Category	PCR microplates
CO Customs tariff code	DE 8479 90 70

Specifications

Ideal fit for*	sciSOURCEPLATE 384 PP BICO PCR-Platte Vollrand, 384 Well, transparent, PP Sarstedt 72.1984, 72.1984.202 PCR-Platte Vollrand, 384 Well, weiß, Low Profile Sarstedt
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	201 g 7.09 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!