

**Product Profile**


Part number	2016-1064
Article name	Adapter for tubes . 96x 0.2 ml
Category	Conical shape / tubes
CO   Customs tariff code	DE   8419 89 98

**Specifications**

Ideal fit for*	
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	210 g   7.41 oz
Recommended max rpm**	2500

\* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

\*\* The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!