

## Adapter specification

---

**Product no.:** 2016-1064

**Product name:** Adapter | 96 0.2ml tubes

---

**Description:** Adapter for cooling and heating of microplates  
Temperature range -10 to 100 °C  
Compatible with all BioShake/ColdPlate/HeatPlate automation modules  
with single point connection

**Ideal fit for:** Eppendorf Safelock 0.2ml tube # 0030124537, Bio-RAD (#MSP-9601)  
and similar disposables

**Material / Color:** Aluminum anodized



**Only using for:** QInstruments Modules with Art. No. 2016-xxx

---

**Manufacturer name:** Quantifoil Instruments GmbH

**Address:** Loebstedter Str. 101 . 07749 Jena . Germany

---



**Drawing:**

