

Adapter specification

Product no.: 2016-1081

Product name: Adapter – 30x 2.0 ml vials, Ø 13.0 mm

Description: Adapter for cooling and heating of microplates
Temperature range -10 to 100 °C
Compatible with all BioShake/ColdPlate/HeatPlate automation modules
with single point connection



Ideal fit for: Cryogenic Vials Corning® #430659, labelled

Material / Color: Aluminum anodized

Only using for: QInstruments Modules with Art. No. 2016-xxx

Manufacturer name: Quantifoil Instruments GmbH

Address: Loebstedter Str. 101 . 07749 Jena . Germany



Drawing:

