

Product Profile



Part number	2016-1121
Article name	Adapter Deep Well Plate . Eppendorf 96/1000 µl
Category	Deep well & storage plates
CO Customs tariff code	DE 8419 89 98

Specifications

Ideal fit for*	<p>Deepwell Plate 96/1000µL Eppendorf 0030503244, 0030 501.209, 0030 501.217, 0030 501.233, 0030 501.241, 0030 505.204, 0030 502.205</p> <p>Abgene 0,8ml storage plate U bottom round PP Thermo AB-0765, AB0765, AB0859, AB-0859</p> <p>Riplate® RW, PS, 1ml Ritter </p>
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	Qinstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	172 g 6.07 oz
Recommended max rpm**	2200

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!