



**Product Profile**

Part number	2016-1121
Article name	Adapter Deep Well Plate . Eppendorf 96/1000 µl
Category	Deep well & storage plates
CO   Customs tariff code	DE   8479 90 70

**Specifications**

Ideal fit for*	<p>Deepwell Plate 96/1000µl Eppendorf   0030503244, 0030 501.209, 0030 501.217, 0030 501.233, 0030 501.241, 0030 505.204, 0030 502.205</p> <p>Abgene 0,8ml storage plate U bottom round PP Thermo   AB-0765, AB0765, AB0859, AB-0859</p> <p>Riplate® RW, PS, 1ml Ritter  </p>
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	172 g   6.07 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!