

Product Profile



Part number	2016-1151
Article name	Adapter- NUNC® & Axygen® deep well 96/2.0 ml
Category	Deep well & storage plates
CO Customs tariff code	DE 8419 89 98

Specifications

Ideal fit for* Material / Color Temperature range Compatible with Mounting material Adapter dimensions Weight Recommended max rpm**	NUNC 96 DW 1ml Thermo 2.0 mL Deep Well Plate for Storage and Collection Corning Nunc 1.0 ml DeepWell™ Plate w shared well technolo Thermo 260251 Aluminum anodized -20 to 100 °C Qinstruments Devices with Art No 2016-xxx (All instruments with single point connection) 1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8 see drawing 141 g 4.97 oz 2200
--	--

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!