

Adapter specification

Product no.: 2016-1151

Product name: Adapter | Deep well, Nunc & Axygen 96/2.0 ml

Description: Adapter for cooling and heating of microplates
Temperature range -10 to 100 °C
Compatible with all BioShake/ColdPlate/HeatPlate automation modules
with single point connection



Ideal fit for: Nunc® 96/2000 µl, #278743, #27875
Axygen® 96/2000 µl, #P-DW-20-C, round bottom

Material / Color: Aluminum anodized


Only using for: QInstruments Modules with Art. No. 2016-xxx

Manufacturer name: Quantifoil Instruments GmbH

Address: Loebstedter Str. 101 . 07749 Jena . Germany



Drawing:



Only for use with QInstruments Modules art. no. 2016-xxxx
Screw 1x BN20146-M3x6-A2 (QI-ID 002255)

