

Adapter specification

Product no.: 2016-1171

Product name: Adapter | Abgene® 96-well storage plate 2.2 ml


Description: Adapter for cooling and heating of microplates
Temperature range -10 to 100 °C
Compatible with all BioShake/ColdPlate/HeatPlate automation modules
with single point connection

Ideal fit for: Abgene® #AB-09032 and similar plates

Material / Color: Aluminum anodized, black matt

Only using for: QInstruments Modules with Art. No. 2016-xxx

Manufacturer name: Quantifoil Instruments GmbH
Address: Loebstedter Str. 101 . 07749 Jena . Germany



Drawing:

