

Product Profile



Part number	2016-1201
Article name	Adapter for microplate . 96 well V-bottom
Category	Standard
CO Customs tariff code	DE 8419 89 98

Specifications

<p>Ideal fit for*</p> <p>Mikrotestplatte, 96 Well, Bodenform: konisch, PS Sarstedt 82.1583, 82.1583.200, 82.1583.100</p> <p>96-well Expanded Volume, PP, V-Bottom Corning </p> <p>Riplate® 96 SRW magnetic, 0.2 ml Ritter </p>	<p>Material / Color</p> <p>Aluminum anodized</p> <p>Temperature range</p> <p>-20 to 100 °C</p> <p>Compatible with</p> <p>QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)</p> <p>Mounting material</p> <p>1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8</p> <p>Adapter dimensions</p> <p>see drawing</p> <p>Weight</p> <p>99 g 3.49 oz</p> <p>Recommended max rpm**</p> <p>2200</p>
--	---

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!