

Adapter specification

Product no.: 2016-1201

Product name: Adapter for microplate . 96 well V-bottom

Description: Adapter for cooling and heating of microplates
Temperature range -10 to 100 °C
Compatible with all BioShake/ColdPlate/HeatPlate automation modules
with single point connection

Ideal fit for:


Material / Color: Aluminum anodized

Only using for: QInstruments Modules with Art. No. 2016-xxxx



Mounting materials: 1x Cylinder screw (BN20146 M3x6-A2, lower head)
1x Torx screwdriver TX8

Manufacturer name: Quantifoil Instruments GmbH

Address: Loebstedter Str. 101 . 07749 Jena . Germany



Drawing:



A-A (1 : 1)

