

Product Profile



Part number	2016-1214
Article name	Adapter - Storage DWP 96 Well 2.0ml
Category	Deep well & storage plates
CO Customs tariff code	DE 8419 89 98

Specifications

<p>Ideal fit for*</p> <p>Material / Color</p> <p>Temperature range</p> <p>Compatible with</p> <p>Mounting material</p> <p>Adapter dimensions</p> <p>Weight</p> <p>Recommended max rpm**</p>	<p>Deep Well MegaBlock®, 96 Well, 1,2 ml, PS, Rundbod Sarstedt 82.1971.002 in PP</p> <p>Riplate® 96 SRW magnetic, 2 ml Ritter </p> <p>Corning® 96-well Clear V-Bottom DW 2 mL PP Corning </p> <p>Aluminum anodized</p> <p>-20 to 100 °C</p> <p>QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)</p> <p>1x Cylinder screw (BN20146 M3x6-A2, lower head) 1x Torx screwdriver TX8</p> <p>see drawing</p> <p>96 g 3.39 oz</p> <p>2200</p>
---	---

* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

** The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!