



**Product Profile**

Part number	2016-1214
Article name	Adapter - Storage DWP 96 Well 2.0ml
Category	Deep well & storage plates
CO   Customs tariff code	DE   8479 90 70

**Specifications**

Ideal fit for*	<p>Deep Well MegaBlock®, 96 Well, 1,2 ml, PS, Rundbod Sarstedt   82.1971.002 in PP</p> <p>Riplate® 96 SRW magnetic, 2 ml Ritter  </p> <p>Corning® 96-well Clear V-Bottom DW 2 mL PP Corning   3960</p>
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	96 g   3.39 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!