

Adapter specification

Product no.: 2016-1215
Product name: Adapter 24 Well Uniplate




Description: Adapter for cooling and heating of microplates
Temperature range -10 to 100 °C
Compatible with all BioShake/ColdPlate/HeatPlate automation modules
with single point connection

Ideal fit for: Whatman™ 7701-5102 UNIPLATE™ 24-Well and similar disposables

Material / Color: Aluminum anodized

Only using for: QInstruments Modules with Art. No. 2016-xxx

Manufacturer name: Quantifoil Instruments GmbH
Address: Loebstedter Str. 101 . 07749 Jena . Germany



Drawing:

