

**Product Profile**



Part number	2016-1217
Article name	Adapter for 96/450µl microplate
Category	Standard
CO   Customs tariff code	DE   8419 89 98

**Specifications**

<p>Ideal fit for*</p> <p>Material / Color</p> <p>Temperature range</p> <p>Compatible with</p> <p>Mounting material</p> <p>Adapter dimensions</p> <p>Weight</p> <p>Recommended max rpm**</p>	<p>Riplate® 96 SRW magnetic, 2 ml Ritter  </p> <p>96 DW PP MASTERBLOCK®, 1 ml Greiner Bio One   780201, 780203, 780204, 780205, 780206, 780215, 780261, 780263, 780264, 780265, 780266</p> <p>96 Square Deep Well Microplate, KingFisher™ Style Azenta  </p> <p>Aluminum anodized</p> <p>-20 to 100 °C</p> <p>Qinstruments Devices with Art No 2016-xxx (All instruments with single point connection)</p> <p>1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8</p> <p>see drawing</p> <p>124 g   4.37 oz</p> <p>2500</p>
---	---

\* Latest information and possibly additional adapters are available online: <https://www.qinstruments.com/service/downloads/>

\*\* The maximal speed depends on a set of parameters like load and height of labware. Please refer to the operation manual of your device for additional information .



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!