

Adapter specification

Product no.: 2016-1217

Product name: Microplate adapter – 96/450µl microplate

Description: Adapter for cooling and heating of microplates
Temperature range -10 to 100 °C
Compatible with all BioShake/ColdPlate/HeatPlate automation modules
with single point connection




Ideal fit for: NUNC® #249944, Storage Plate, PP, 450µl, round and similar disposables

Material / Color: Aluminum anodized

Only using for: QInstruments Modules with Art. No. 2016-xxx

Manufacturer name: Quantifoil Instruments GmbH

Address: Loebstedter Str. 101 . 07749 Jena . Germany



Drawing:

