

Adapter specification

Product no.: 2016-1218

Product name: Adapter for Apricot® 96/2.2ml Deep Well Plate

Description: Adapter for cooling and heating of microplates
Temperature range -10 to 100 °C
Compatible with all BioShake/ColdPlate/HeatPlate automation modules
with single point connection



Ideal fit for: Apricot DWP-22-96-SQ-U-C and similar disposables

Material / Color: Aluminum anodized

Only using for: QInstruments modules with art. no. 2016-xxx

Manufacturer name: Quantifoil Instruments GmbH

Address: Loebstedter Str. 101 . 07749 Jena . Germany



Drawing:

