


**Product Profile**

Part number	2016-1219
Article name	Adapter - 24 Deep Well Plate
Category	Conical shape / tubes
CO   Customs tariff code	DE   8479 90 70

**Specifications**

Ideal fit for*	KingFisher™ Plastics for 24 deep-well format Thermo   95040480, 97002580, 95040470, 97002610, 97002620, 95040470B, 97002610B
Material / Color	Aluminum anodized
Temperature range	-20 to 100 °C
Compatible with	QInstruments Devices with Art No 2016-xxx (All instruments with single point connection)
Mounting material	1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8
Adapter dimensions	see drawing
Weight	177 g   6.24 oz
Recommended max rpm	The maximum speed depends on a set of parameters like used device, liquid volume, height of labware, ... A general applicable value cannot be specified. Please refer to the operation manual of your device for additional information.

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!