


**Product Profile**

|                          |                                   |
|--------------------------|-----------------------------------|
| Part number              | 2016-1315                         |
| Article name             | Adapter - 5x Tecan 25ml Reservoir |
| Category                 | Reservoir                         |
| CO   Customs tariff code | DE   8479 90 70                   |

**Specifications**

|                     |  |
|---------------------|--|
| Ideal fit for*      |  |
| Material / Color    | Aluminum anodized  |
| Temperature range   | -20 to 100 °C  |
| Compatible with     | ColdPlate, ColdPlate slim and HeatPlate, instruments with mixing function are not supported. |
| Mounting material   | 1x Cylinder screw (BN20146 M3x6-A2, lower head)   1x Torx screwdriver TX8                    |
| Adapter dimensions  | see drawing  |
| Weight              | 331 g   11.68 oz   |
| Recommended max rpm | 0 (instruments with mixing function are not supported)                                       |

\* Latest information and possibly additional adapters are available on request: <https://www.qinstruments.com/contact/>



TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE!